

**Protocol MSU\_MSMC\_007**, Version 1.1

## **Methoximation and trimethylsilylation of amino acids and organic acids for GC/MS analysis**

Last updated on Nov 17, 2020

### **Reagents**

Methoxyamine hydrochloride (Sigma-Aldrich #89803)

Pyridine (dry), (Sigma-Aldrich #270970-4X25ML)

*N*-Methyl-*N*-trimethylsilyltrifluoroacetamide (MSTFA) containing 1% trimethylsilyl chloride (TMSCl), in sealed glass ampules (Sigma-Aldrich #69478)

### **Supplies**

1.7-mL polypropylene microcentrifuge tubes with locking caps (VWR #490016-245)

Calibrated 1000- $\mu$ L pipetter and pipet tips

Calibrated 10- $\mu$ L pipetter and pipet tips

Vortexer

Ultrasound water bath

Labeled amber autosampler vials (BMB Stores #21140 with low volume (250  $\mu$ L) glass inserts (BMB Stores #51832085) and PTFE-lined screw caps (BMB Stores #06718904)

Oven or heated block

Analytical balance (to 0.1 mg)

Spatula, precleaned

### **Samples**

Use extracts of one of the following (after evaporation of solvents to dryness, typically in a screw cap vial; can use autosampler vials with inserts; extraction details are in separate SOPs):

blood serum or plasma (30  $\mu$ L)

urine (50-100  $\mu$ L)

Cell cultures ( $10^7$  cells)

Cell culture media (50  $\mu$ L)

Homogenized tissue (plant, animal, or other): 2-25 mg of tissue

### **Procedure**

1. Set oven (or heated block) temperature to 60°C.
2. Prepare labels for microcentrifuge tubes and GC vials
3. Use a spatula to weigh 0.040 g of methoxyamine hydrochloride into a 1.7 mL microcentrifuge tube labeled as "40 mg/mL methoxyamine-HCl in pyridine"
4. Transfer 1000  $\mu$ L of dry pyridine into the above microcentrifuge tube; seal the tube.
5. Vortex briefly, then ultrasonicate for 15 minutes; ensure that all of the solid has dissolved before proceeding.
6. Add 100  $\mu$ L of the methoxyamine-HCl/pyridine solution to each dried sample, blank, QC sample, and calibration standard.
7. Heat tubes at 60°C for 12-24 hours; allow tubes to cool to room temperature.
8. Add 100  $\mu$ L of MSTFA + 1% TMSCl and seal the tube.
9. Heat tubes at 60°C for 12-24 hours; allow tubes to cool to room temperature.
10. Transfer 50  $\mu$ L of each reaction mixture to an autosampler vial equipped with low-volume insert; cap the vial and transfer it to the GC/MS autosampler.
11. Transfer the remaining reaction materials to a separate labeled autosampler vial (no insert necessary) for storage.