

Protocol MSU_MSMC_007a

Gas chromatography-mass spectrometry profiling of central metabolites derivatized by methoximation/trimethylsilylation (GC/MS operating procedures)

Materials

Gas chromatography-mass spectrometry instrument (e.g. Agilent 7890 GC coupled to Agilent 5975 (or 7010) mass spectrometer with autosampler (e.g. CTC PAL)
Retention index mix of n-alkanes, C10-C40, all even, 50 ng/μL (Sigma-Aldrich #68281-2ML-F)
Sample metadata should be provided with each sample

Experiment design

Each batch of samples to be analyzed should include the following:

One vial containing retention index mix of n-alkanes (50 ng/μL)

One vial (for each 20 samples) containing a procedural blank prepared with derivatizing agents (methoxyamine hydrochloride in pyridine and MSTFA + 1% TMSCI (see Protocol MSU_MSMC_007)

Metabolite extracts prepared and derivatized according to Protocol MSU_MSMC_007 and randomized in order of analysis (e.g. use RAND function in Microsoft Excel)

One vial containing hexanes (solvent blank)

Order of sample analyses

1. Hexanes solvent blank
2. Procedural blank
3. n-alkanes retention mix standard solution
4. Derivatized metabolite extracts

Experimental parameters for GC/MS analysis

Column:	Agilent VF-5ms + 10 m EZ-Guard (30 m long, 0.25 mm I.D., 0.25 μm film; 10 m uncoated retention gap); part #CP9013
Mobile phase:	Helium
Mobile phase flow rate	1.2 mL/min (flow controlled)
Injection volume	1.0 μL
Injector type	Mixed mode injector (MMI) or programmed temperature vaporizer (PTV)
Injector temperature (program)	50°C ramped to 250°C at 12°/min; hold at 250° until 50 minutes post-injection
Injection type	Splitless; purge at 0.4 min; purge flow 5 mL/min
Column temperature	50°C (1 min hold), increase at 6°/min to 330°; hold at 320° for 5 minutes (total acquisition time: 52.67 minutes)
Ionization method	Electron ionization (70 eV)
Ion source temperature	250°C
Transfer line temperature	280°C

Data acquisition type	Scan
Scan range	<i>m/z</i> 80-600
Solvent delay (min)	3.0
Scan time (ms)	250
Detector gain	1
Mass calibration procedure	Tune must have been evaluated within the past seven days
<i>Autosampler procedures</i>	
PAL Cycle	GC-Liq4-V3
Wash solvent 1	Ethyl acetate
Wash solvent 2	Hexanes
Pre-clean with solvent 1	3
Pre-clean with solvent 2	3
Preclean with sample	1
% Syringe fill for cleaning	80
Syringe fill volume	8 μ L
Sample amount for cleaning (μ L)	1
Filling speed (μ L/s)	5
Filling strokes	6
Volume for filling strokes	6
Injection speed (μ L/s)	50
Pre-inject delay (ms)	500
Post-inject delay (ms)	500
Post clean with solvent 1	3
Post clean with solvent 2	3

Data processing to be described in a separate protocol.