

Protocol MSU_MSMC_014

Procedure for preparation of nucleotide phosphate extracts for UHPLC/MS/MS quantification using ion-pairing separation on a reversed phase column

(following the protocol of Fernandez and Waters, *Current Protocols in Microbiology*, 52, e74. doi: 10.1002/cpmc.74)

Last modified August 31, 2020

Materials

Microbial culture at OD₆₀₀ = 0.5.

Speed Vac with cold trap.

-20°C freezer.

Refrigerated centrifuge with rotor for 1.5-mL microcentrifuge tubes.

15-mL Falcon polypropylene centrifuge tubes.

Crushed ice.

N,N-Dimethylhexylamine (Sigma-Aldrich cat# 308102-5G).

c-di-GMP fluorinated (InVivoGen STING Ligand - Difluoro 3'3'-c-di-GMP, cyclic [FdG(3',5')pFdG(3',5')p]; Catalog #tlrl-nacdgf); abbreviated below as c-di-GMPF.

Acetic acid (glacial).

Methanol, HPLC grade.

Acetonitrile, HPLC grade.

Water, MilliQ (> 18 MΩ).

Formic acid (98%) (Fisher Scientific cat# 60-006-17).

LB media.

Preparation of reagents and solvents

1. Nucleotide extraction solution (0.1 M formic acid in acetonitrile/methanol/water 40:40:20 v/v/v) containing 25 nM c-di-GMPF as internal standard:
 - 1) Add 4000 μL of HPLC grade acetonitrile to a 15-mL polypropylene centrifuge tube.
 - 2) Add 4000 μL of HPLC grade methanol to the 15-mL polypropylene centrifuge tube.
 - 3) Add 2000 μL of MilliQ water to the 15-mL polypropylene centrifuge tube.
 - 4) Label acetonitrile/methanol/water (40:40:20 v/v/v) as solvent A.
 - 5) Mix 25 μL of 10 μM c-di-GMPF and 38.5 μL of 98% formic acid and 9936.5 μL of solvent A in a 15-mL polypropylene centrifuge tube.
 - 6) Label the tube as “Nucleotide Extraction Solution with 25 nM IS” and store it at -20°C.
2. Mobile phase solvent A:
 - 1) Fill 1000 mL of MilliQ water in a 1-liter HPLC solvent bottle.
 - 2) Add 1.387 mL of N,N-Dimethylhexylamine (DMHA).
 - 3) Add 164.7 μL of acetic acid; these ingredients give 8 mM DMHA, 2.8 mM acetic acid.
 - 4) Cap the solvent bottle and place it in an ultrasound bath. Ultrasonicate for 15 minutes. Allow solution to stand overnight before use.

- 5) Label the solvent bottle (8 mM DMHA/2.8 mM acetic acid in water, pH~9).

Preparing internal standard c-di-GMPF at 50 nM:

1. Take 250 μL of 1 μM stock IS c-di-GMPF, add 4500 μL of MP-A, to make 50 nM of IS solution.
2. Store it at -20°C (or -80°C) until ready to use by LC/MS/MS

Preparing standard curve standard solutions:

1. From 1 μM standard stock solution, use MP-A to make serial dilutions to obtain concentrations (nM): 500, 250, 100, 50, 20, 10, 5, 2.
2. Then dilute each standard stock solution with an equal volume of 50 nM internal standard c-di-GMPF solution.
3. The final concentration of internal standard c-di-GMPF is then 25 nM.

500 nM —>	250 nM
250 nM —>	125 nM
100 nM —>	50 nM
50 nM —>	25 nM
20 nM —>	10 nM
10 nM —>	5 nM
5 nM —>	2.5 nM
2 nM —>	1 nM

Sample preparation procedure

1. Chill centrifuge to 4°C .
2. Prepare 100 μL of extraction solution per sample, store on ice or in a -20°C explosion-proof or flammable storage freezer.
3. Transfer two replicate 1000- μL aliquots of $\text{OD}_{600} = 0.5$ culture into separate 1.5-mL polypropylene microcentrifuge tubes.
4. Centrifuge tubes at 15,000 $\times g$ for 30 seconds.
5. During the above centrifugation step, transfer 500 μL from the original culture to a new tube and add 500 μL of LB medium; measure and record OD_{600}
 - a. This OD_{600} measurement is to be used for normalizing the nucleotide concentrations to the number of cells extracted
 - b. Alternatively, concentrations can be normalized to protein content of extracted sample
 - c. More detailed information on normalization can be found in the Fernandez and Waters protocol.
6. From the original tubes (from step 4), quickly remove the supernatant using a 1000- μL pipette.

7. Resuspend the pellet in 100 μ L of cold nucleotide extraction solution (with 25 nM IS) to lyse cells and quench metabolism.
8. Place the quenched pellet/extraction solution tubes in a -20°C freezer for 20 minutes.
9. Remove tubes from freezer and centrifuge for 15 minutes at 15,000 x g, 4°C .
10. Transfer supernatant to a new 1.5-mL microcentrifuge tube.
11. Evaporate extracts to dryness under vacuum (Speed Vac) without application of heat.
12. Redissolve in 100 μ L of mobile phase solvent A.
13. Transfer to HPLC autosampler vials with 200- μ L glass insert.

Note: Redissolving samples in mobile phase A containing N,N-Dimethylhexylamine (MDHA) is critical for optimal ion-pairing and reduces the chances of split peaks during chromatography.